

Labycare-Autopsy Saw With Vacuum Dust Collector (SKSS-AP017)



Autopsy Saw With Vacuum Dust Collector Manufacturer by LABYCARE

Description :

An Autopsy Saw with Vacuum Dust Collector is designed for safe, precise bone cutting during post-mortem examinations. The integrated suction system captures bone dust and aerosols, reducing contamination and improving operator safety in forensic and pathology environments.

Specifications :

Type: Electric oscillating autopsy saw

Dust Control: Built-in vacuum dust collector

Blade: Stainless steel, interchangeable

Speed: Variable control

Handle: Ergonomic, anti-slip

Power Supply: 220–240 V AC

Application: Mortuaries, forensic and pathology labs

labycare

www.labycare.com | Email: labycare@gmail.com