

Labycare-Anthropometric Rods (SKSS-AP018)



Anthropometric Rods Manufacturer by LABYCARE

Description :

Anthropometric Rods are precision instruments used to measure body height, limb length, and skeletal dimensions during anthropological, forensic, and medical examinations. Designed for accuracy and durability, they are essential for post-mortem analysis, research, and educational use.

Specifications :

Material: Aluminum / Stainless steel

Measurement Range: Up to 200 cm

Scale: Clear, engraved metric markings

Sections: Telescopic / fixed segments

Accuracy: ± 1 mm

Design: Lightweight, portable

Application: Forensic labs, anthropology, mortuaries

labycare

www.labycare.com | Email: labycare@gmail.com