

Labycare-Bullet Probe Instruments (SKSS-AP019)



Bullet Probe Instruments Manufacturer by LABYACRE

Description :

Bullet Probe Instruments are used during forensic and post-mortem examinations to trace bullet paths and locate foreign bodies within tissues. They provide accurate trajectory assessment while minimizing tissue damage, aiding medico-legal investigations.

Specifications :

Material: Stainless steel

Design: Smooth, rounded probe tip

Length: Various sizes available

Finish: Corrosion-resistant, polished

Sterilization: Autoclavable

Application: Forensic pathology, autopsy rooms

labycare

www.labycare.com | Email: labycare@gmail.com