

Labycare-Headrest Plastic Six Position (SKSS-AP020)



Headrest Plastic Six Position Manufacturer by LABYCARE

Description :

Plastic Six Position Headrest is designed to support and stabilize the head during autopsy and post-mortem procedures. Adjustable height and angle ensure proper positioning, comfort, and safety while maintaining hygiene and ease of cleaning in mortuary environments.

Specifications :

Material: High-grade medical plastic

Positions: Six adjustable settings

Design: Anatomical, contoured shape

Surface: Smooth, easy to clean

Load Capacity: Suitable for adult use

Application: Autopsy tables, mortuaries, forensic labs

