

Labycare-Mail Gloves SS Chain (SKSS-AP021)



Mail Gloves SS Chain Manufacturer by LABYCARE

Description :

Stainless Steel Chain Mail Gloves provide superior hand protection against cuts, punctures, and sharp instruments during autopsy, meat processing, and industrial tasks. Ergonomic design ensures flexibility, safety, and long-lasting performance in high-risk working environments.

Specifications :

Material: Stainless steel chain mail

Design: Interlinked rings, flexible fit

Cuff: Adjustable strap

Size: Small / Medium / Large

Protection: Cut and puncture resistant

Hygiene: Easy to clean, corrosion resistant

Application: Autopsy rooms, mortuaries, food processing units

labycare

www.labycare.com | Email: labycare@gmail.com