

Labycare-Osteometric Boards (SKSS-AP022)



Osteometric Boards Manufacturer by LABYCARE

Description :

Osteometric Boards are precision measuring tools used in anthropology and forensic science to determine bone length and skeletal dimensions accurately. Designed for reliability and ease of use, they are essential for stature estimation and skeletal analysis.

Specifications :

Material: Hardwood / Aluminum / Acrylic

Measurement Range: Up to 60–100 cm

Scale: Clear metric graduations

Design: Fixed vertical end and sliding horizontal end

Accuracy: ± 1 mm

Surface: Smooth, easy to clean

Application: Forensic labs, anthropology departments, medical colleges

labycare

www.labycare.com | Email: labycare@gmail.com