

Labycare-Silicone Oil (SKSS-AP024)



Silicone Oil Manufacturer by LABYCARE

Description :

Silicone Oil is a clear, inert, and high-performance fluid widely used for lubrication, heat transfer, damping, and laboratory applications. It offers excellent thermal stability, low volatility, and water repellency, ensuring reliable performance across a wide temperature range.

Specifications :

Chemical Type: Polydimethylsiloxane (PDMS)

Appearance: Clear, colorless liquid

Viscosity: Available in multiple grades

Temperature Range: -50°C to $+250^{\circ}\text{C}$

Properties: Non-toxic, non-reactive, water repellent

Packaging: Bottle / can / drum

Application: Laboratory, pharmaceutical, industrial use

labycare

www.labycare.com | Email: labycare@gmail.com