

# Labycare-Mortuary Cabinet 1 Compartment (SKSSMC1)



Mortuary Cabinet 1 Compartment Manufacturer by **LABYCARE**

## **Description :**

Mortuary Chamber are made of high grade Stainless Steel 304 Grade

Mortuary Chamber is also known as Morque Cabinet, Mortuary Chamber and Morque Chamber

These Cabinets comprises of one compartment with carriage for storage of cadaver.

Compartment is equipped with a hinged door and locking arrangement.

The Insulating layer is made of imported puf with high thermal retardation against cold loss

The high quality ultra-mute totally enclosed type has individual compressors for each body

The refrigeration system of each chamber is individually controlled and temperature can be adjusted without interference

The High quality gel strip coated door gasket with strong aging resistance ability lengthens the service life of the door gasket

## **Specifications :**

Usage/Application	Hospital , Autopsy
Material	SS304
Brand	LABYCARE
Color	Silver
Dimension	Internal size: 1950
Weight	150 kg

Temperature Range Degree Celsius	2 deg to 8 deg
Warranty	1 year
Is It Lockable	yes
Order Quantity	1
Voltage	220V
Country of Origin	Made in India

#### Technical Specifications:

- **Temperature Range:** The One Body Mortuary Chamber maintains a temperature range of 2 to 8 degrees Celsius, ensuring optimal preservation conditions.
- **Construction:** Built with high-quality materials, including stainless steel interior walls and floors, the Mortuary Chamber is durable and easy to clean.
- **Insulation:** The Mortuary Chamber is insulated with high-density polyurethane foam to ensure efficient temperature control and energy conservation.
- **Cooling System:** The Mortuary Chamber is equipped with a reliable and energy-efficient cooling system, ensuring consistent temperature levels.
- **Power Requirements:** The Mortuary Chamber operates on standard electrical power, with voltage and frequency requirements varying based on the specific model and configuration.
- **Surface Finish:** The Mortuary Chamber's surface finish meets international standards for hygiene and durability, ensuring a clean and professional appearance.
- **Doors:** The Mortuary Chamber is equipped with a single door, designed for easy access and equipped with safety features such as a door alarm and emergency release mechanism.
- **Interior Layout:** The Mortuary Chamber features stainless steel trays or racks for the storage of bodies, providing efficient use of space and easy organization.
- **Controls:** Temperature and other settings are easily controlled via a user-friendly interface, allowing for precise adjustment to meet specific requirements.
- **Safety Features:** In addition to a door alarm and emergency release mechanism, the Mortuary Chamber may include other safety features such as a temperature alarm and power failure alarm.
- **Customization:** The Mortuary Chamber is available in various sizes to accommodate the needs of different facilities, with options for additional features and accessories based on specific requirements.

---

labycare

[www.labycare.com](http://www.labycare.com) | Email: [labycare@gmail.com](mailto:labycare@gmail.com)