

Labycare-Hydraulic Dead Body Lifter (SKSS-AP031)



Hydraulic Dead Body Lifter Manufacturer by LABYCARE

Description :

The hydraulic dead body lifter is designed for safe, smooth, and dignified handling of deceased bodies. Featuring a sturdy stainless steel frame and foot-operated hydraulic system, it allows effortless height adjustment, reduced manual strain, and stable transfer in hospitals, mortuaries, and forensic facilities.

Specifications :

Operation: Foot-operated hydraulic pump

Material: Stainless steel / powder-coated steel frame

Height adjustment: Smooth, variable hydraulic lift

Load capacity: Up to 200 kg (approx.)

Castors: Heavy-duty swivel castors with brakes

Top: Removable stainless steel tray/platform

Finish: Rust-resistant, easy to clean

Application: Mortuary and hospital use
