

Labycare-Electro Hydraulic Dead Body Lifter (SKSS-AP032)



Electro Hydraulic Dead Body Lifter Manufacturer by LABYCARE

Description :

The electro hydraulic dead body lifter is designed for effortless and dignified handling of deceased bodies. Equipped with an electric motor and hydraulic system, it enables smooth height adjustment at the push of a button, ensuring safety, stability, and efficiency in mortuaries and hospitals.

Specifications :

- Operation: Electric motor with hydraulic lifting system
- Control: Push-button or hand-held remote
- Material: Stainless steel / powder-coated steel frame
- Height adjustment: Smooth, variable electro-hydraulic lift
- Load capacity: Up to 200–250 kg (approx.)
- Castors: Heavy-duty swivel castors with brakes
- Power supply: 220–240 V AC, 50 Hz

- Finish: Corrosion-resistant, easy to clean

labycare

www.labycare.com | Email: labycare@gmail.com