

Labycare-Blood Bag Tube Sealer Without Battery Backup (SKSS-BBP03)



Blood Bag Tube Sealer Without Battery Backup Manufacturer by LABYCARE

Description :

A Blood Bag Tube Sealer Without Battery Backup is a device used to create secure, leak-proof seals on blood bag tubing during collection or processing. It operates on mains power, offering reliable performance for safe handling, storage, and transport of blood products in hospitals, blood banks, and clinical laboratories.

Specifications :

- Mains-powered operation (no battery backup)
- Produces secure, leak-proof seals on tubing
- Compact, ergonomic design
- Easy to operate with minimal training
- Durable construction for frequent use
- Ideal for blood banks, hospitals, and labs

labycare

www.labycare.com | Email: labycare@gmail.com