

Labycare-Blood Bag Triple CPD (SKSS-BBP029)



Blood Bag Triple CPD Manufacturer by LABYCARE

Description :

A Blood Bag Triple CPD is a sterile, single-use blood collection system with three interconnected bags containing CPD anticoagulant solution. It allows collection and separation of whole blood into components like plasma and red cells while preserving quality and viability for safe storage and transfusion in blood banks and hospitals.

Specifications :

- Triple-bag system with CPD anticoagulant
- Sterile, medical-grade single-use design
- Leak-proof tubing with sealed ports
- Facilitates efficient component separation
- Compatible with blood processing equipment
- Suitable for blood banks and clinical use

labycare

www.labycare.com | Email: labycare@gmail.com