

Labycare-Transfer Bag (SKSS-BBP034)



Transfer Bag Manufacturer by LABYCARE

Description :

A Transfer Bag is a sterile, single-use medical bag used to transfer blood components or fluids from one container to another during processing. It ensures safe, contamination-free handling, supports component separation, and maintains product integrity in blood banks, hospitals, and clinical laboratories.

Specifications :

- Sterile, single-use medical-grade bag
- Made from high-strength, flexible plastic
- Leak-proof tubing with secure ports
- Compatible with blood bag systems
- Available in various capacities
- Suitable for blood banks and clinical use

labycare

www.labycare.com | Email: labycare@gmail.com