

Labycare-Plasma Extractor Electric (SKSS-BBP036)



Plasma Extractor Electric Manufacturer by LABYCARE

Description :

An Electric Plasma Extractor is an automated device used in blood banks to separate plasma from centrifuged whole blood. It applies controlled pressure using an electric motor, ensuring uniform, efficient plasma expression with minimal cell damage, improved accuracy, and enhanced workflow compared to manual extraction methods.

Specifications :

- Electric motor-driven plasma extraction
- Controlled and uniform pressure application
- Compatible with standard blood bag systems
- Easy-to-operate control panel
- Durable, easy-to-clean construction
- Suitable for blood banks and clinical laboratories

labycare

www.labycare.com | Email: labycare@gmail.com