

Labycare-Vaccine Carriers 2.90 ltr (SKSS-CCE011)



Vaccine Carriers 2.90 ltr Manufacturer by LABYCARE

Description :

A Vaccine Carrier 2.90 Ltr is an insulated container designed to safely transport vaccines while maintaining the required cold chain temperature. Used with ice packs, it preserves vaccine potency during immunization sessions, outreach programs, and short-distance medical transport in healthcare and public health settings.

Specifications :

- 2.90-liter storage capacity
- High-quality thermal insulation for temperature retention
- Compatible with standard ice packs
- Lightweight and easy to carry
- Durable, easy-to-clean outer body
- Suitable for immunization and vaccination programs

labycare

www.labycare.com | Email: labycare@gmail.com