

Labycare-Dental Impression Tray Standard (SKSS-DP021)



Dental Impression Tray Standard Manufacturer by LABYCARE

Description :

A Dental Impression Tray Standard is used to hold and support impression material while taking accurate dental impressions of teeth and oral tissues. It ensures proper shape, stability, and uniform material distribution, helping dentists obtain precise models for diagnosis, treatment planning, and fabrication of dental prosthetics.

Specifications :

Made from high-quality stainless steel or medical-grade plastic

Standard design suitable for routine dental impressions

Smooth edges for patient comfort

Rigid construction for accurate impression results

Easy to clean and disinfect / autoclavable (metal type)

Available in multiple sizes for upper and lower arches

labycare

www.labycare.com | Email: labycare@gmail.com