

Labycare-Height measuring stand (SKSSHMS919)



Height measuring stand

Description :

Height measuring stand

Specifications :	Height measuring stand	Height Measuring Stand : Telescopic stand, adjustable built on dial type weighing scale.
TECHNICAL SPECIFICATION:		• Height measuring scale, floor model without weighing scale. Measuring range: 20 - 210cm with 1mm graduation. Measuring in inches and centimeters. Height stand disassembles into 4 sections and easy to transport in disassembled condition • Suitable for height measuring at any location. Accurate, portable, durable and light weight with all plastic parts • Suitable for medical practices, hospitals, nursing home, school, home. Supplied with knock down condition with how to assemble manual in english language.

