

## Labycare-Automated Tissue Processor (SKSSATP921)



### Automated Tissue Processor

#### **Description :**

Automated Tissue Processor

|                                |                            |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
|--------------------------------|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b><u>Specifications :</u></b> | Automated Tissue Processor | Histokennette Automatic Tissue Processor with safety device. Equipped with latest micro switch mechanism. Supplied with one stainless steel basket rotor, one stainless steel basket, 24 tissue capsules, 10 one litre beaker, 2 thermostatically controlled wax baths, 6 calibrated time discs, one cutting pliers, dust cover and instruction manual. The safety device checks whether the wax is melted or not and hence protects the tissues by not allowing them to shift from 10th to 11th or 12th stations if the wax is not melted. |
|                                | Brand                      | SKSS                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |

---

labycare

www.labycare.com | Email: labycare@gmail.com