


Labycare-Room Thermometer Tabletop (SKSS-GMP072)



Room Thermometer Tabletop Manufacturer by LABYCARE

Description :

Displays the temperature, time and humidity all at the same time.

Clock and calendar function (month and date).

12-hour or 24-hour display

Large Character LCD display

Temperature range: -10~ + 50 °C (14~122 °F)

Humidity range: 10%~99%RH

We are Exporters and Suppliers of Room Thermometer Tabletop at wholesale competitive prices from India.

Specifications :

Displays the temperature, time and humidity all at the same time.

Clock and calendar function (month and date).

12-hour or 24-hour display

Large Character LCD display

Temperature range: -10~ + 50 °C (14~122 °F)

Humidity range: 10%~99%RH

We are Exporters and Suppliers of Room Thermometer Tabletop at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com