

Labycare-Mayo Tray Stainless Steel 202 Grade (SKSS-HH020)



Mayo Tray Stainless Steel 202 Grade Manufacturer by LABYCARE

Description :

Mayo Trays are commonly known as Mayo Table Tray, Mayo Instrument Tray and Mayo Surgical Tray.

Mayo Tray is manufactured with stainless steel 202 grade without having any sharp edges.

These trays have seamless jointless construction.

Mayo Trays are supplied without cover

Mayo Tray is used above the Mayo Trolley for keeping Surgical Instruments

Stainless Steel Mayo Tray is durable and Reusable.

We are Manufacturer, Exporters and Suppliers of Mayo Tray Stainless Steel 202 Grade at wholesale competitive prices from India.

Specifications :

Mayo Trays are commonly known as Mayo Table Tray, Mayo Instrument Tray and Mayo Surgical Tray.

Mayo Tray is manufactured with stainless steel 202 grade without having any sharp edges.

These trays have seamless jointless construction.

Mayo Trays are supplied without cover

Mayo Tray is used above the Mayo Trolley for keeping Surgical Instruments

Stainless Steel Mayo Tray is durable and Reusable.

We are Manufacturer, Exporters and Suppliers of Mayo Tray Stainless Steel 202 Grade at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com