

Labycare-Instruments Box Stainless Steel 304 Grade (SKSS-HH082)



Instruments Box Stainless Steel 304 Grade Manufacturer by LABYCARE

Description :

Instrument Box is also known as Surgical Box

It is used to store Surgical Instruments.

Instrument Box is made of Stainless steel 304 Grade

These boxes are highly polished Stainless steel and come with friction fit lid

The instruments box makes it easier for surgeons, to take surgical instruments during operations quickly. This saves their time during a serious operation.

Our Instrument boxes are high quality and durable.

We are Manufacturer, Exporters and Suppliers of Instruments Box Stainless Steel 304 Grade at wholesale competitive prices from India.

Specifications :

Instrument Box is also known as Surgical Box

It is used to store Surgical Instruments.

Instrument Box is made of Stainless steel 304 Grade

These boxes are highly polished Stainless steel and come with friction fit lid

The instruments box makes it easier for surgeons, to take surgical instruments during operations quickly. This saves their time during a serious operation.

Our Instrument boxes are high quality and durable.

We are Manufacturer, Exporters and Suppliers of Instruments Box Stainless Steel 304 Grade at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com