

Labycare-Scrub Station Single Tap (SKSS-HH096)



Scrub Station Single Tap Manufacturer by LABYCARE

Description :

Scrub Station Single is also known as Single Bay Scrub Station and Scrub Sink Single.

Single Bay Scrub Station is attached with infrared photo sensor.

Digital timer has defined timer interval 1,3,5 and 10 minutes.

1 No. Mixture valve and 2 Nos. Solenoid valves are provided along with manual option foot / knees operated mechanism .

The Anti-splash metal body construction is made of stainless steel 304 Grade with polished finished thickness 1.5 mm.

The Scrub Station Single Tap is provided with front access panel for ease of maintenance and repairs

We are Manufacturer, Exporters and Suppliers of Scrub Station Single Tap at wholesale competitive prices from India.

Specifications :

Scrub Station Single is also known as Single Bay Scrub Station and Scrub Sink Single.

Single Bay Scrub Station is attached with infrared photo sensor.

Digital timer has defined timer interval 1,3,5 and 10 minutes.

1 No. Mixture valve and 2 Nos. Solenoid valves are provided along with manual option foot / knees

operated mechanism .

The Anti-splash metal body construction is made of stainless steel 304 Grade with polished finished thickness 1.5 mm.

The Scrub Station Single Tap is provided with front access panel for ease of maintenance and repairs

We are Manufacturer, Exporters and Suppliers of Scrub Station Single Tap at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com