

Labycare-Scrub Station Wall Mounted Double Tap (SKSS-HH0100)



Scrub Station Wall Mounted Double Tap Manufacturer by LABYCARE

Description :

Scrub Station Wall Mounted Double is also known as Double Bay Scrub Station Wall Mounted and Scrub Station Sink Wall Mounted Double.

Double Bay Scrub Station Wall Mounted is attached with infrared photo sensor.

Digital timer has defined timer interval 1,3,5 and 10 minutes.

2 No. Mixture valve is provided along with manual Elbow operated mechanism .

The Anti-splash metal body construction is made of stainless steel 304 Grade with polished finished thickness 1.5 mm.

We are Manufacturer, Exporters and Suppliers of Scrub Station Double Tap at wholesale competitive prices from India.

Specifications :

Scrub Station Wall Mounted Double is also known as Double Bay Scrub Station Wall Mounted and Scrub Station Sink Wall Mounted Double.

Double Bay Scrub Station Wall Mounted is attached with infrared photo sensor.

Digital timer has defined timer interval 1,3,5 and 10 minutes.

2 No. Mixture valve is provided along with manual Elbow operated mechanism .

The Anti-splash metal body construction is made of stainless steel 304 Grade with polished finished thickness 1.5 mm.

We are Manufacturer, Exporters and Suppliers of Scrub Station Double Tap at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com