

Labycare-Scrub Station Wall Mounted Triple Tap (SKSS-HH0101)



Scrub Station Wall Mounted Triple Tap Manufacturer by LABYCARE

Description :

Scrub Station Wall Mounted Triple is also known as Triple Bay Scrub Station Wall Mounted and Scrub Station Sink Wall Mounted Triple.

Triple Bay Scrub Station Wall Mounted is attached with infrared photo sensor.

Digital timer has defined timer interval 1,3,5 and 10 minutes.

3 No. Mixture valve is provided along with manual Elbow operated mechanism .

The Anti-splash metal body construction is made of stainless steel 304 Grade with polished finished thickness 1.5 mm.

We are Manufacturer, Exporters and Suppliers of Scrub Station Triple Tap at wholesale competitive prices from India.

Specifications :

Scrub Station Wall Mounted Triple is also known as Triple Bay Scrub Station Wall Mounted and Scrub Station Sink Wall Mounted Triple.

Triple Bay Scrub Station Wall Mounted is attached with infrared photo sensor.

Digital timer has defined timer interval 1,3,5 and 10 minutes.

3 No. Mixture valve is provided along with manual Elbow operated mechanism .

The Anti-splash metal body construction is made of stainless steel 304 Grade with polished finished thickness 1.5 mm.

We are Manufacturer, Exporters and Suppliers of Scrub Station Triple Tap at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com