

Labycare-Rotarod (SKSSR929)



Rotarod 6 Compartments for Rat and Mice

Description :

Rotarod 6 Compartments for Rat and Mice

<u>Specifications :</u>	Rota Rod Compartments	6	Digital Rota Rod apparatus incorporating electronic stepless variable speeds from 3 RPM to 25 RPM of a common serrated 50 mm. shaft. Non Transparent (black) perspex compartments walls with independent electronic digital timers and contactless photoelectric time recording operation supplied with interchangeable 25mm. dia. Rod for Mice and 50 mm. dia rod for Rat - 6 compartments.
	Brand		SKSS

labycare

www.labycare.com | Email: labycare@gmail.com