

Labycare-X-Ray Thyroid Shield Collar (SKSS-XRP011)



X-Ray Thyroid Shield Collar Manufacturer by LABYCARE

Description :

A thyroid guard or thyroid shield is a lead collar that wraps around your neck to block the radiation that's generated in making X-ray images. During an X-ray, the majority of radiation needed to create the images goes exactly where it's aimed. In the case of a mammogram, most radiation exposure occurs in your breast.

X Ray Thyroid Shield Collar is finest collar design that fits thyroid gland perfectly

Thyroid Shield Collar has wide coverage

Thyroid Shield Collar comes in one size that can fit all

Thyroid Shield Collar is light weight which helps in easy movement of neck.

We are Exporters and Suppliers of X-Ray Thyroid Shield Collar at wholesale competitive prices from India

Specifications :

A thyroid guard or thyroid shield is a lead collar that wraps around your neck to block the radiation that's generated in making X-ray images. During an X-ray, the majority of radiation needed to create the images goes exactly where it's aimed. In the case of a mammogram, most radiation exposure occurs in your breast.

X Ray Thyroid Shield Collar is finest collar design that fits thyroid gland perfectly

Thyroid Shield Collar has wide coverage

Thyroid Shield Collar comes in one size that can fit all

Thyroid Shield Collar is light weight which helps in easy movement of neck.

We are Exporters and Suppliers of X-Ray Thyroid Shield Collar at wholesale competitive prices from India

labycare

www.labycare.com | Email: labycare@gmail.com