

Labycare-X-Ray Lead Cap (SKSS-XRP012)



X-Ray Lead Cap Manufacturer by LABYCARE

Description :

The X-Ray Cap is a reusable cap that is designed to protect the wearer's head from scatter radiation. It is made from durable fabrics to ensure maximum levels of comfort during procedures, and features an elasticated band to provide a good degree of security.

X-ray lead caps provide x-ray protection for the top of the head.

Lead equivalent is .5mm for the sides and .25mm for the top.

Lead Caps are fitted with a Velcro closure.

Size and diameter is measured with tape around head just above the ears.

We are Exporters and Suppliers of X-Ray Lead Cap at wholesale competitive prices from India.

Specifications :

The X-Ray Cap is a reusable cap that is designed to protect the wearer's head from scatter radiation. It is made from durable fabrics to ensure maximum levels of comfort during procedures, and features an elasticated band to provide a good degree of security.

X-ray lead caps provide x-ray protection for the top of the head.

Lead equivalent is .5mm for the sides and .25mm for the top.

Lead Caps are fitted with a Velcro closure.

Size and diameter is measured with tape around head just above the ears.

We are Exporters and Suppliers of X-Ray Lead Cap at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com