

Labycare-X-Ray Face Shield Mask (SKSS-XRP023)



X-Ray Face Shield Mask Manufacturer by LABYCARE

Description :

X-Ray Face mask with full visor, shaped to protect the face.

The two side knobs make it possible to adjust the mask to any angle in order to reduce exposure to scattered radiation that could filter up from below. There is also a knob at the back to regulate the circumference to ensure a better fit

Face Shield Mask frame material is made of Plastic

Lead equivalence Front Shield - 0.10mm Lead Acrylic

Face Shield Mask is regular frame size with adjustable head pieces.

We are Exporters and Suppliers of X-Ray Face Shield Mask at wholesale competitive prices from India.

Specifications :

X-Ray Face mask with full visor, shaped to protect the face.

The two side knobs make it possible to adjust the mask to any angle in order to reduce exposure to scattered radiation that could filter up from below. There is also a knob at the back to regulate the circumference to ensure a better fit

Face Shield Mask frame material is made of Plastic

Lead equivalence Front Shield - 0.10mm Lead Acrylic

Face Shield Mask is regular frame size with adjustable head pieces.

We are Exporters and Suppliers of X-Ray Face Shield Mask at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com