

Labycare-X-Ray Cassette Pass Box (SKSS-XRP025)



X-Ray Cassette Pass Box Manufacturer by LABYCARE

Description :

X-Ray Cassette Pass Box is made of heavy Gauge Mild Steel with double wall construction

All joints are welded smoothly

X-Ray Cassette Pass Box has separate compartment labeled "Exposed" and "Unexposed"

X-Ray Cassette Pass Box is Rugged, Reliable and Durable

X-Ray Cassette Pass Box has Automatic Lock System

X-Ray Cassette Pass Box comes with Lead Lining

X-Ray Cassette Pass Box is used for Medical X-Ray Film, Ultrasound Film, Dry View Film.

We are Manufacturer, Exporters and Suppliers of X-Ray Cassette Pass Box at wholesale competitive prices from India.

Specifications :

X-Ray Cassette Pass Box is made of heavy Gauge Mild Steel with double wall construction

All joints are welded smoothly

X-Ray Cassette Pass Box has separate compartment labeled "Exposed" and "Unexposed"

X-Ray Cassette Pass Box is Rugged, Reliable and Durable

X-Ray Cassette Pass Box has Automatic Lock System

X-Ray Cassette Pass Box comes with Lead Lining

X-Ray Cassette Pass Box is used for Medical X-Ray Film, Ultrasound Film, Dry View Film.

We are Manufacturer, Exporters and Suppliers of X-Ray Cassette Pass Box at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com