

Labycare-X-Ray Cassette Film Cart (SKSS-XRP026)



X-Ray Cassette Film Cart Manufacturer by LABYCARE

Description :

Cassette Film Cart is made of mild steel with baked enamel finish.

Cassette/Film Cart has six 6" wide compartments on the top shelf vertical clearance.

Lower shelf is modified to accept 14"x17" cassettes.

Rubber mat isn fitted inside each compartment.

Side Loading Carts enhance file room efficiency and organization with flexible applications.

Sturdy 4" (102mm) plate castors: two stationary, two swivel.

Padded bar-style handle and scuff resistant top.

We are Manufacturer, Exporters and Suppliers of X-Ray Cassette Film Cart at wholesale competitive prices from India.

Specifications :

Cassette Film Cart is made of mild steel with baked enamel finish.

Cassette/Film Cart has six 6" wide compartments on the top shelf vertical clearance.

Lower shelf is modified to accept 14"x17" cassettes.

Rubber mat isn fitted inside each compartment.

Side Loading Carts enhance file room efficiency and organization with flexible applications.

Sturdy 4" (102mm) plate castors: two stationary, two swivel.

Padded bar-style handle and scuff resistant top.

We are Manufacturer, Exporters and Suppliers of X-Ray Cassette Film Cart at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com