

Labycare-X-Ray Protection Suspension Ceiling (SKSS-XRP032)



X-Ray Protection Suspension Ceiling Manufacturer by LABYCARE

Description :

Protection Suspension Ceiling Type comes with Lead Rubber

The transparent part of the shield is made of Lead Acrylic Material

Shield Size : 400mm (H) x 600mm (W)

Curtain Size : 350mm (H) x 600mm (W)

Lead Equivalent is 0.50mmPb

We are Manufacturer, Exporters and Suppliers of X-Ray Protection Suspension Ceiling at wholesale competitive prices from India.

Specifications :

Protection Suspension Ceiling Type comes with Lead Rubber

The transparent part of the shield is made of Lead Acrylic Material

Shield Size : 400mm (H) x 600mm (W)

Curtain Size : 350mm (H) x 600mm (W)

Lead Equivalent is 0.50mmPb

We are Manufacturer, Exporters and Suppliers of X-Ray Protection Suspension Ceiling at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com