

Labycare-X-Ray Ceiling Suspension Shield with O.T. Light (SKSS-XRP034)



X-Ray Ceiling Suspension Shield with O.T. Light Manufacturer by LABYCARE

Description :

Illuminance (Lux) - 60,000~1,20,000

Color Temperature -4000+ 500

Spot Diameter (mm) -160~280

Illuminance Depth(mm) -1000

Brightness Adjustment - 8-Stage continuous light adjustment

Total Irradiance - 580w/m²

The ratio between irradiance

And illuminance - 407mW/m²lx

Color reduction index(Ra)—93

Power Supply Voltage - ~220v/50hz

Input Power (W) - 200

24V Bulb Power (W) - 150

The transparent part of the shield is made of Lead Acrylic Material

Shield Size : 750mm (H) x 600mm (W)

Lead Equivalent is 0.50mmPb

We are Manufacturer, Exporters and Suppliers of X-Ray Ceiling Suspension Shield with O.T. Light at wholesale competitive prices from India.

Specifications :

Illuminance (Lux) - 60,000~1,20,000

Color Temperature -4000+ 500

Spot Diameter (mm) -160~280

Illuminance Depth(mm) -1000

Brightness Adjustment - 8-Stage continuous light adjustment

Total Irradiance - 580w/m²

The ratio between irradiance

And illuminance - 407mW/m²lx

Color reduction index(Ra)—93

Power Supply Voltage - ~220v/50hz

Input Power (W) - 200

24V Bulb Power (W) - 150

The transparent part of the shield is made of Lead Acrylic Material

Shield Size : 750mm (H) x 600mm (W)

Lead Equivalent is 0.50mmPb

We are Manufacturer, Exporters and Suppliers of X-Ray Ceiling Suspension Shield with O.T. Light at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com