

## **Labycare-X-Ray View Box Single with Sensor & Dimming Function (SKSS-XRP049)**



X-Ray View Box Single with Sensor & Dimming Function Manufacturer by LABYCARE

### **Description :**

X-Ray View Box is also known as X-Ray Film Viewer

X-Ray View box's body is made of Mild Steel Cabinet with Frame and with Acrylic Sheet

X-Ray View Box is for Single X-Ray Film

X-Ray View Box is powered by 12 V, 50-60Hz adapter.

X-Ray View Box is Sensor based, Edge Light will ON once you insert a film and it will OFF after removing the same.

X-Ray View box has heavy duty switch & socket with Shock Proof Mechanism

X-Ray View Box has plastic clips to hold the X-ray film.

X-Ray View Box is light weight & portable

Dimmer Switch is provided to adjust the intensity of Light.

We are Manufacturer, Exporters and Suppliers of X-Ray View Box Single with Sensor & Dimming Function at wholesale competitive prices from India.

### **Specifications :**

X-Ray View Box is also known as X-Ray Film Viewer

X-Ray View box's body is made of Mild Steel Cabinet with Frame and with Acrylic Sheet

X-Ray View Box is for Single X-Ray Film

X-Ray View Box is powered by 12 V, 50-60Hz adapter.

X-Ray View Box is Sensor based, Edge Light will ON once you insert a film and it will OFF after removing the same.

X-Ray View box has heavy duty switch & socket with Shock Proof Mechanism

X-Ray View Box has plastic clips to hold the X-ray film.

X-Ray View Box is light weight & portable

Dimmer Switch is provided to adjust the intensity of Light.

We are Manufacturer, Exporters and Suppliers of X-Ray View Box Single with Sensor & Dimming Function at wholesale competitive prices from India.

---

labycare

[www.labycare.com](http://www.labycare.com) | Email: [labycare@gmail.com](mailto:labycare@gmail.com)