

## Labycare-X-Ray Cassette Holder Chest Stand Wall Mounted (SKSS-XRP064)



X-Ray Cassette Holder Chest Stand Wall Mounted Manufacturer by LABYCARE

### **Description :**

X Ray stand supports in stabilizing the mammography system, leading to the accurate positioning of X Ray screen, thus exact recording of X Ray images.

Chest Stand allows to accommodate the screen plates and cassettes of varied sizes and gets strong hold of the devices.

X-Ray Cassette Holder Chest Stand Wall Mounted is easily foldable

X-Ray Cassette Holder Chest Stand Wall Mounted is Adjustable and

Corrosion resistance

Chest Stand is easily mounted to the wall with 4 screws

The Cassette holding channels slide freely along the vertical rails yet lock firmly at any elevation by means of two positive spring locks/ manual nobs

The cassette holder takes cassettes size upto 14"x17"

All parts are durable metal, finished in powder coating except the rods

We are Manufacturer, Exporters and Suppliers of X-Ray Cassette Holder Chest Stand Wall Mounted at wholesale competitive prices from India.

### **Specifications :**

X Ray stand supports in stabilizing the mammography system, leading to the accurate positioning of X Ray screen, thus exact recording of X Ray images.

Chest Stand allows to accommodate the screen plates and cassettes of varied sizes and gets strong hold of the devices.

X-Ray Cassette Holder Chest Stand Wall Mounted is easily foldable

X-Ray Cassette Holder Chest Stand Wall Mounted is Adjustable and

Corrosion resistance

Chest Stand is easily mounted to the wall with 4 screws

The Cassette holding channels slide freely along the vertical rails yet lock firmly at any elevation by means of two positive spring locks/ manual nobs

The cassette holder takes cassettes size upto 14"x17"

All parts are durable metal, finished in powder coating except the rods

We are Manufacturer, Exporters and Suppliers of X-Ray Cassette Holder Chest Stand Wall Mounted at wholesale competitive prices from India.

---

labycare

[www.labycare.com](http://www.labycare.com) | Email: [labycare@gmail.com](mailto:labycare@gmail.com)