

Labycare-Colorimeter (SKSSC935)



Digital Photo Colorimeter 8 Filters mains operated with optical density in 2½ Digit LED Display, 1ml solution measurement.

Description :

Digital Photo Colorimeter 8 Filters mains operated with optical density in 2½ Digit LED Display, 1ml solution measurement.

Specifications : TECHNICAL

SPECIFICATION:

Display	2½ Digit	2½ Digit
Range	400 nm to 700 nm	400 nm to 700 nm
Filters	Glass Filters	Glass Filters
Accuracy	± 2% of F.S.	± 2% of F.S.
Light Sources	6V . 10 w	6V . 10 w
Sensor	Selenium Photocell/Diode	Selenium Photocell/Diode
Power Supply	220 VAC ± 10%	220 VAC ± 10%
Battery	-----	-----
Filters	8 Filters	5 Filters

labycare

www.labycare.com | Email: labycare@gmail.com