

Labycare-Digital Blood Cell Counter or Digital Differential Counter (SKSS-PLP04)



Digital Blood Cell Counter or Digital Differential Counter _____

Description :

Differential Blood Cell Counter has microprocessor controlled system which offers accurate and quick differential counts.

Digital Blood Cell Counter is featured with reliability, accuracy, stability and is a sleek device. These cell counters are available with 11 operating keys which are suitable for clinical labs for showing the primary bacterial or virus infections.

Cell Counting Key - 8 Key

Total Key - 11 Key (8 Cell keys and 3 Function keys)

Body - Metal Body

Display - 16 X 2 Line LCD Display

Main Input - 230V \pm 10% AC, 50 Hz

Dimension - 250 x 230 x 150 mm (LxBxH)

Weight - 1.5 kg (Approx.)

1.16 Digits Alphanumeric LCD Display

Automatic percentage (%) calculation

Individual cells counting facility

In-built Stop Watch and Timer function

Audible beep facility

We are Manufacturer, Exporters and Suppliers of Digital Blood Cell Counter or Digital Differential Counter at wholesale competitive prices from India.

Specifications :

Differential Blood Cell Counter has microprocessor controlled system which offers accurate and quick differential counts.

Digital Blood Cell Counter is featured with reliability, accuracy, stability and is a sleek device. These cell counters are available with 11 operating keys which are suitable for clinical labs for showing the primary bacterial or virus infections.

Cell Counting Key - 8 Key

Total Key - 11 Key (8 Cell keys and 3 Function keys)

Body - Metal Body

Display - 16 X 2 Line LCD Display

Main Input - 230V \pm 10% AC, 50 Hz

Dimension - 250 x 230 x 150 mm (LxBxH)

Weight - 1.5 kg (Approx.)

1.16 Digits Alphanumeric LCD Display

Automatic percentage (%) calculation

Individual cells counting facility

In-built Stop Watch and Timer function

Audible beep facility

We are Manufacturer, Exporters and Suppliers of Digital Blood Cell Counter or Digital Differential Counter at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com