

Labycare-Pipette Filler 3 Way Rubber (SKSS-MRP028)



Pipette Filler 3 Way Rubber Manufacturer by LABYCARE

Description :

Pipette Filler 3 Way is made of superior natural rubber

Pipette Filler has three valve safety pipette filler

It eliminates the need for mouth pipetting.

These Pipette Fillers are chemical resistant

It is easy to use and has one hand operation. Automatic valve is removable for easy cleaning

Pipet filler can be used for all pipets from 0.05 mL to 100 mL .

Simply squeeze the ball and it is ready to use. Press to fill, press to empty, then press to eject the pipet.

We are Manufacturer, Exporters and Suppliers of Pipette Filler 3 Way Rubber at wholesale competitive prices from India.

Specifications :

Pipette Filler 3 Way is made of superior natural rubber

Pipette Filler has three valve safety pipette filler

It eliminates the need for mouth pipetting.

These Pipette Fillers are chemical resistant

It is easy to use and has one hand operation. Automatic valve is removable for easy cleaning

Pipet filler can be used for all pipets from 0.05 mL to 100 mL .

Simply squeeze the ball and it is ready to use. Press to fill, press to empty, then press to eject the pipet.

We are Manufacturer, Exporters and Suppliers of Pipette Filler 3 Way Rubber at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com