

Labycare-Pipette Filler 3 Way Silicone (SKSS-MRP029)



Pipette Filler 3 Way Silicone Manufacturer by LABYCARE

Description :

Pipette Filler 3 Way is made of superior silicone

Pipette Filler has three valve safety pipette filler

It eliminates the need for mouth pipetting.

Silicone provides greater chemical resilience than natural rubber.

Silicone has greater chemical resistance than natural rubber.

It is easy to use and has one hand operation. Automatic valve is removable for easy cleaning

Pipet filler can be used for all pipets from 0.05 mL to 100 mL .

Simply squeeze the ball and it is ready to use. Press to fill, press to empty, then press to eject the pipet.

We are Manufacturer, Exporters and Suppliers of Pipette Filler 3 Way Silicone at wholesale competitive prices from India.

Specifications :

Pipette Filler 3 Way is made of superior silicone

Pipette Filler has three valve safety pipette filler

It eliminates the need for mouth pipetting.

Silicone provides greater chemical resilience than natural rubber.

Silicone has greater chemical resistance than natural rubber.

It is easy to use and has one hand operation. Automatic valve is removable for easy cleaning

Pipet filler can be used for all pipets from 0.05 mL to 100 mL .

Simply squeeze the ball and it is ready to use. Press to fill, press to empty, then press to eject the pipet.

We are Manufacturer, Exporters and Suppliers of Pipette Filler 3 Way Silicone at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com