

Labycare-Anaerobic apparatus (SKSSAA887)



Anaerobic apparatus

Description :

Anaerobic apparatus

Specifications : Anaerobic apparatus	Anaerobic Culture Jar made of stainless steel. Jar size 23cm. X 15 cm. With top flange recessed for a toro dial rubber sealing ring and a needle valve for evacuation and a sachet containing a catalyst. A side arm on the jar connecting to the glass tube containing semi solid Lucas indicator. The jar is supplied without semi solid Lucas indicator. With pressure gauge.
Brand	SKSS