

Labycare-Endotracheal Tubes Plain PVC (SKSS-MD01)



Endotracheal Tubes Plain PVC Manufacturer by LABYCARE

Description :

Endotracheal Tubes is also known as ET Tube.

Endotracheal Tube or ET tube has a smooth inside/outside surface for ease in intubation, extubation and suctioning.

Endotracheal Tube or ET tube has smooth polished murphy eye, distal end and tapered balloon end for atraumatic intubation and extubation.

ET tube is kink resistant inflation tube to ensure patient safety during cuff inflation and deflation.

Radio opaque line and markings on the ET tube facilitate identification of tube position.

Endotracheal tube has universal 15 mm connector at the proximal end.

ET tube is E.T.O Sterile and pyrogen free.

It is meant for single patient use only.

We are Manufacturer, Exporters and Suppliers of Endotracheal Tubes Plain PVC at wholesale competitive prices from India.

Specifications :

Endotracheal Tubes is also known as ET Tube.

Endotracheal Tube or ET tube has a smooth inside/outside surface for ease in intubation, extubation and suctioning.

Endotracheal Tube or ET tube has smooth polished murphy eye, distal end and tapered balloon end

for atraumatic intubation and extubation.

ET tube is kink resistant inflation tube to ensure patient safety during cuff inflation and deflation.

Radio opaque line and markings on the ET tube facilitate identification of tube position.

Endotracheal tube has universal 15 mm connector at the proximal end.

ET tube is E.T.O Sterile and pyrogen free.

It is meant for single patient use only.

We are Manufacturer, Exporters and Suppliers of Endotracheal Tubes Plain PVC at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com