

Labycare-Endotracheal Tube Plain South Oral Preformed (SKSS-MD03)



Endotracheal Tube Plain South Oral Preformed Manufacturer by LABYCARE

Description :

Endotracheal Tube is also known as ET Tube.

ET Tube is used to maintain a patent airway and mechanically ventilate patients.

Endotracheal Tube is made from non-toxic, non-irritant medical grade PVC.

Endotracheal Tube has uniform spiral metal reinforced tube wall to prevent the tube from kinking.

Reinforcement supports tube position and verification during X-ray.

Accurate and safe positioning ensured by graduation marks and tip-to-tip radiopaque line.

UET Tube is fitted with universal 15 mm connector at proximal end.

Recommended use of bite block / Guedel airway.

E.T.O sterile and pyrogen free.

ET Tube is for single patient use.

We are Manufacturer, Exporters and Suppliers of Endotracheal Tube Plain South Oral Preformed at wholesale competitive prices from India.

Specifications :

Endotracheal Tube is also known as ET Tube.

ET Tube is used to maintain a patent airway and mechanically ventilate patients.

Endotracheal Tube is made from non-toxic, non-irritant medical grade PVC.

Endotracheal Tube has uniform spiral metal reinforced tube wall to prevent the tube from kinking.

Reinforcement supports tube position and verification during X-ray.

Accurate and safe positioning ensured by graduation marks and tip-to-tip radiopaque line.

UET Tube is fitted with universal 15 mm connector at proximal end.

Recommended use of bite block / Guedel airway.

E.T.O sterile and pyrogen free.

ET Tube is for single patient use.

We are Manufacturer, Exporters and Suppliers of Endotracheal Tube Plain South Oral Preformed at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com