

Labycare-OT Light Single Dome Ceiling 21 LED (SKSS-OL01)



OT Light Single Dome Ceiling 21 LED Manufacturer by LABYCARE

Description :

A Single Dome Ceiling-Mounted OT Light with 21 LEDs provides bright, shadow-free illumination for minor and general surgical procedures, ensuring clear visibility and color accuracy in the operating theatre.

Specifications :

- **Type:** Ceiling-mounted, single dome
- **Light Source:** 21 high-intensity LEDs
- **Illumination:** Approx. 80,000–120,000 Lux
- **Color Temperature:** 4,000–4,500 K (natural white)
- **Light Field Diameter:** 160–200 mm
- **Shadow Reduction:** Multi-faceted reflector / LED arrangement
- **Control:** Intensity control via handle or wall panel
- **Heat Emission:** Minimal, cool light
- **Power Supply:** 220–240 V AC
- **Applications:** Minor OT, general surgery, procedure rooms

labycare

www.labycare.com | Email: labycare@gmail.com