

Labycare-OT Light Single Dome Ceiling 30 LED (SKSS-OL02)



OT Light Single Dome Ceiling 30 LED Manufacturer by LABYCARE

Description :

A Single Dome Ceiling-Mounted OT Light with 30 LEDs delivers powerful, uniform, and shadow-free illumination for general and minor surgical procedures, ensuring excellent visibility and true color rendering in the operating theatre.

Specifications :

Type: Ceiling-mounted, single dome

Light Source: 30 high-intensity LEDs

Illumination: Approx. 1,00,000–1,60,000 Lux

Color Temperature: 4,000–4,500 K (natural white)

Light Field Diameter: 160–220 mm

Shadow Reduction: Advanced multi-LED optical design

Control: Intensity adjustment via control panel / sterile handle

Heat Emission: Cool light, low IR radiation

Power Supply: 220–240 V AC, 50/60 Hz

Applications: General surgery, minor OT, procedure rooms

labycare

www.labycare.com | Email: labycare@gmail.com