

Labycare-OT Light Double Dome Ceiling 60 LED (SKSS-OL07)



OT Light Double Dome Ceiling 60 LED Manufacturer by LABYCARE

Description :

A double dome ceiling OT light with 60 LEDs provides bright, shadow-free illumination for major surgical procedures. It ensures uniform light distribution, excellent color rendering, and adjustable intensity, making it ideal for operation theatres, emergency rooms, and surgical suites.

Specifications :

- Light source: 60 high-intensity LEDs (30 + 30)
- Illumination: Up to 1,60,000 lux (each dome)
- Color temperature: 4,000–4,500 K
- Mounting: Ceiling mounted, double dome
- Control: Adjustable intensity & focus
- Power supply: 220–240 V AC

labycare

www.labycare.com | Email: labycare@gmail.com