

Labycare-Autopsy Tables (SKSSAT378)



Autopsy Tables

Description :

Autopsy Tables

<u>Specifications :</u> Autopsy Tables	AUTOPSY TABEL :-All S.S. Construction sturdy support framework. One endmounted intergrated large SS wash basin with one free rotating faucet and one long reach hand shower. Built in under table plumbing and drainage. (External pluming not included).Comes with one SS tubular portable and adjustabel overhead light on casters.
Brand	SKSS

labycare

www.labycare.com | Email: labycare@gmail.com