

## Labycare-OT Light Double Dome Ceiling 104 LED (SKSS-OL09)



OT Light Double Dome Ceiling 104 LED Manufacturer by LABYCARE

### **Description :**

A double dome ceiling OT light with 104 LEDs provides high-intensity, shadowless illumination for advanced surgical procedures. It ensures uniform light distribution, excellent depth of illumination, and accurate color rendering, making it suitable for major and super-specialty operation theatres.

### **Specifications :**

- Light source: 104 high-intensity LEDs (52 + 52)
- Illumination: Up to 2,00,000 lux (each dome)
- Color temperature: 4,000–4,500 K
- Mounting: Ceiling mounted, double dome
- Control: Digital control panel with adjustable intensity
- Power supply: 220–240 V AC

---

labycare

[www.labycare.com](http://www.labycare.com) | Email: [labycare@gmail.com](mailto:labycare@gmail.com)