

## Labycare-OT Light Mobile 21 LED (SKSS-OL011)



OT Light Mobile 21 LED Manufacturer by LABYCARE

### **Description :**

A Mobile OT Light with 21 LEDs is designed for minor surgeries and examination procedures, offering bright, uniform, and shadow-free illumination with easy mobility and flexible positioning.

### **Specifications :**

**Type:** Mobile / Floor-mounted

**Light Source:** 21 high-intensity LEDs

**Illumination:** Approx. 60,000–100,000 Lux

**Color Temperature:** 4,000–4,500 K (natural white)

**Light Field Diameter:** 150–200 mm

**Arm Movement:** Flexible arm with smooth positioning

**Control:** Intensity control on light head / control panel

**Mobility:** Heavy-duty base with lockable castors

**Heat Emission:** Cool light, minimal IR & UV radiation

**Power Supply:** 220–240 V AC, 50/60 Hz

**Applications:** Minor OT, emergency rooms, clinics, procedure rooms

---

labycare

[www.labycare.com](http://www.labycare.com) | Email: [labycare@gmail.com](mailto:labycare@gmail.com)