

Labycare-OT Light Mobile 30 LED (SKSS-OL012)



OT Light Mobile 30 LED Manufacturer by LABYCARE

Description :

A Mobile OT Light with 30 LEDs provides powerful, uniform, and shadow-free illumination for minor to general surgical procedures, with flexible positioning and easy mobility across operation theatres.

Specifications :

- **Type:** Mobile / Floor-mounted
- **Light Source:** 30 high-intensity LEDs
- **Illumination:** Approx. 80,000–1,20,000 Lux
- **Color Temperature:** 4,000–4,500 K (natural white)
- **Light Field Diameter:** 160–220 mm
- **Arm Movement:** Flexible, spring-balanced arm
- **Control:** Intensity adjustment on light head / control panel
- **Mobility:** Stable base with antistatic lockable castors
- **Heat Emission:** Cool light, minimal IR & UV radiation
- **Power Supply:** 220–240 V AC, 50/60 Hz

- **Applications:** Minor OT, emergency, clinics, procedure rooms

labycare

www.labycare.com | Email: labycare@gmail.com