

Labycare-OT Light Mobile 48 LED (SKSS-OL013)



OT Light Mobile 48 LED Manufacturer by LABYCARE

Description :

A Mobile OT Light with 48 LEDs delivers high-intensity, uniform, and shadow-free illumination for general and major surgical procedures, with smooth maneuverability and precise light positioning.

Specifications :

Type: Mobile / Floor-mounted

Light Source: 48 high-intensity LEDs

Illumination: Approx. 1,20,000–1,80,000 Lux

Color Temperature: 4,000–4,500 K (natural white)

Light Field Diameter: 180–230 mm

Color Rendering Index (CRI): ?95

Arm Movement: Spring-balanced, flexible arm

Control: Intensity control on light head / control panel

Mobility: Heavy-duty base with antistatic lockable castors

Heat Emission: Cool light, minimal IR & UV radiation

Power Supply: 220–240 V AC, 50/60 Hz

Applications: General OT, emergency OT, trauma, specialty procedures

labycare

www.labycare.com | Email: labycare@gmail.com