

Labycare-OT Light Mobile 52 LED (SKSS-OL014)



OT Light Mobile 52 LED Manufacturer by LABYCARE

Description :

A Mobile OT Light with 52 LEDs provides ultra-bright, uniform, and shadow-free illumination for general and major surgical procedures, offering flexibility, precise positioning, and easy mobility.

Specifications :

Type: Mobile / Floor-mounted

Light Source: 52 high-intensity LEDs

Illumination: Approx. 1,40,000–2,00,000 Lux

Color Temperature: 4,000–4,500 K (natural white)

Light Field Diameter: 180–240 mm

Color Rendering Index (CRI): ?95

Arm Movement: Spring-balanced, flexible arm

Control: Intensity control on light head / control panel

Mobility: Stable base with antistatic lockable castors

Heat Emission: Cool light with minimal IR & UV radiation

Power Supply: 220–240 V AC, 50/60 Hz

Applications: Major OT, emergency, trauma & specialty surgeries

labycare

www.labycare.com | Email: labycare@gmail.com